


File : C:\msdchem\1\data\072616 SCAQMD\IOM Mix 1st Run.D  
 Operator : Dave  
 Acquired : 26 Jul 2016 12:59 using AcqMethod SCAQMD M311.M  
 Instrument : GCMS TDU  
 Sample Name : IOM Mix 1st Run  
 Misc Info :  
 Vial Number: 13

GC-MS



C:\msdchem\1\methods\SCAQMD M311.M  
 Wed Jul 27 09:20:01 2016

GC-MS

Control Information

Sample Inlet : GC  
 Injection Source : Manual  
 Mass Spectrometer : Enabled

Injection Location: Front

No Sample Prep method has been assigned to this method.

Oven  
 Equilibration Time 1 min  
 Max Temperature 280 degrees C  
 Slow Fan Disabled  
 Oven Program On  
 40 °C for 2 min  
 then 10 °C/min to 100 °C for 3 min  
 then 10 °C/min to 250 °C for 8 min  
 Run Time 34 min

Sample Overlap  
 Sample overlap is not enabled

Front SS Inlet He  
 Mode Split  
 Heater On 250 °C  
 Pressure On 1.8804 psi  
 Total Flow On 42 mL/min  
 Septum Purge Flow On 3 mL/min  
 Gas Saver On 20 mL/min After 2 min  
 Split Ratio 25 :1  
 Split Flow 37.5 mL/min

Back PTV Inlet He  
 Mode Split  
 Heater Off  
 Pressure Off  
 Total Flow Off  
 Septum Purge Flow Off  
 Gas Saver Off  
 Split Ratio 20 :1  
 Split Flow 0 mL/min  
 Cryo Off

Thermal Aux 2 (MSD Transfer Line)  
 Heater On  
 Temperature Program On  
 250 °C for 0 min  
 Run Time 34 min

Column #1  
 Agilent Restek 13870 Restek Rxi-624SIL MS  
 280 °C: 30 m x 320 µm x 1.8 µm  
 In: Front SS Inlet He  
 Out: Vacuum

|                      |               |
|----------------------|---------------|
| (Initial)            | 40 °C         |
| Pressure             | 1.8804 psi    |
| Flow                 | 1.5 mL/min    |
| Average Velocity     | 44.411 cm/sec |
| Holdup Time          | 1.1258 min    |
| Flow Program         | Off           |
| 1.5 mL/min for 0 min |               |
| Run Time             | 34 min        |

|           |          |
|-----------|----------|
| Signals   |          |
| Test Plot | Save Off |
|           | 50 Hz    |
| Test Plot | Save Off |
|           | 50 Hz    |
| Test Plot | Save Off |
|           | 50 Hz    |
| Test Plot | Save Off |
|           | 50 Hz    |

GERSTEL MAESTRO

SYSTEM SETTINGS

|                   |             |
|-------------------|-------------|
| Maestro Runtime   | : 34.00 min |
| GC Cool Down Time | : 2.00 min  |
| Cryo Timeout      | : 15.00 min |

GERSTEL CIS

CIS : not used - use these parameters if it becomes 'used'

CRYO COOLING

Cryo Cooling : not used

TEMPERATURE PROGRAM

|                     |              |
|---------------------|--------------|
| Heater Mode         | : Standard   |
| Initial Temperature | : 60 °C      |
| Equilibration Time  | : 1.00 min   |
| Initial Time        | : 1.00 min   |
| Ramp 1              |              |
| Rate                | : 10.00 °C/s |
| End Temp            | : 240 °C     |
| Hold Time           | : 1.00 min   |
| Ramp 2              |              |
| Rate                | : 0.0 °C/s   |

GERSTEL MPS Liquid Injection

Syringe : 10ul

SAMPLE PARAMETERS

|                  |              |
|------------------|--------------|
| Sandwich         | : not used   |
| Inj. Volume      | : 1.0 uL     |
| Air Volume below | : 0.0 uL     |
| Inj. Speed       | : 50.00 uL/s |
| Fill Volume      | : 5.0 uL     |
| Fill Strokes     | : 4          |
| Fill Speed       | : 2.50 uL/s  |

Viscosity Delay : 2 s  
Eject Speed : 50.00 uL/s  
Pre Inj. Delay : 0 s  
Post Inj. Delay : 0 s  
Inj. Penetration : 40.00 mm  
Sample Tray Type : VT98  
Vial Penetration : 31.00 mm

CLEANING PARAMETERS

Preclean Sample : 0  
Wash Station 1 : Wash1  
Preclean Solv.1 : 4  
Postclean Solv.1 : 4  
Fill Speed Solv.1 : 5.00 uL/s  
Viscosity Delay Solv.1 : 2 s  
Eject Speed Solv.1 : 50.00 uL/s  
Wash Station 2 : Wash2  
Preclean Solv.2 : 0  
Postclean Solv.2 : 0  
Fill Speed Solv.2 : 5.00 uL/s  
Viscosity Delay Solv.2 : 2 s  
Eject Speed Solv.2 : 50.00 uL/s

MS ACQUISITION PARAMETERS

General Information  
-----

Tune File : atune.u  
Acquisition Mode : Scan

MS Information  
--

Solvent Delay : 0.00 min  
EMV Mode : Relative  
Relative Voltage : 0  
Resulting EM Voltage : 1659

[Scan Parameters]

Low Mass : 20.0  
High Mass : 400.0  
Threshold : 100  
Sample # : 3 A/D Samples 8  
Plot 2 low mass : 30.0  
Plot 2 high mass : 300.0

[MSZones]

MS Source : 230 C maximum 250 C  
MS Quad : 150 C maximum 200 C

Timed Events  
-----

[Timed MS Detector Table Entries]

| Time (min) | State (MS On/Off) |
|------------|-------------------|
| 4.40       | Off               |
| 5.20       | On                |

END OF MS ACQUISITION PARAMETERS