




1

SCHEDULE	
9:00 AM	INTRODUCTIONS
9:00 - 9:45 AM	SURVEY OF M313
9:45 AM- LUNCH	DEMONSTRATION [LAB]
1:00 - 2:30 PM	TELECONFERENCE [CC-2]
2:30 - 4:30 PM	PROCESSING/DISCUSSION [CC-2]
4:30 PM	WHOLE LAB TOUR [LAB]

2



3

INSTRUMENT CONFIGURATION

4

CONSUMABLES SEPTA, LINER, SYRINGE

SPLIT VS. SPLITLESS CHROMATOGRAPHY VS
EQUILIBRATION TIME

POST-COLUMN SPLIT SPLITTER INSTALLATION

INSTRUMENT CONFIGURATION

5

IOM PREPARATION PREP DEMONSTRATION

INTERPRETING RESULTS

INSTRUMENT
RECONFIGURATION LOOP

INSTRUMENT OPTIMIZATION MIX

6

IOM Component	Evaluation Of
C6 - C15	Distribution of Molecular Weights to FID
BFB	Mass Spec Ionization [TO-15]
Triglyme	Instrument Sensitivity
EGDE (IS)	Injector Reliability
Ethylene Glycol	Separation from Internal Standard
Propylene Glycol	
Methyl Palmitate	Endpoint Retention Time

INSTRUMENT OPTIMIZATION MIX

7

IOM Component	QC Requirements
C6 - C15	Mass-adjusted area counts within 85 - 115% (normalized to C10)
BFB	Passes TO-15 Criteria
Triglyme	80 - 120% recovery at 0.1 g/L
EGDE (IS)	50 - 150% recovery
Ethylene Glycol	90% separation of each component (EG/EGDE/PG)
Propylene Glycol	
Methyl Palmitate	Within 0.1 min prior to and following sample analysis

IOM WORKSHEET

8


Pre-Sample IOM												
Container	Standard	CAS	Index #	Manufacturer	Lot #	Density	% purity	Retention time	Area Count	AreaMass	% of Normalized	
Hexane	11543	1	Fisher	13817D 2108257	0.676	99.9	4.340	31496417	0.0652	0.0645	402776529	99.48
Isodecane	14265	3	Sigma-Aldrich	05471R	0.804	99.9	7.615	33008173	0.0652	0.0675	502489194	91.87
Nonane	111842	2	Sigma-Aldrich	ST80788DV	0.724	99.0	14.951	37723294	0.0712	0.0795	135171541	97.86
Decane	124185	2	Sigma-Aldrich	85137H	0.730	99.6	17.32	39687307	0.0729	0.0726	546338146	100.00
Undecane	1120214	3	Sigma-Aldrich	W89029TV	0.740	99.4	19.231	40327146	0.0743	0.0739	546027116	100.45
Dodecane	112403	2	Sigma-Aldrich	W89750AV	0.750	99.0	20.864	41990521	0.0750	0.0742	564388583	103.27
Tridecane	62955	1	Aldrich	W89156RV	0.756	99.0	22.368	43266887	0.0761	0.0761	563848742	103.20
Tetradecane	62954	3	Sigma-Aldrich	ST80371V	0.764	99.0	23.789	43176348	0.0760	0.0762	573843094	100.04
Pentadecane	62953	1	Sigma-Aldrich	101170J	0.769	99.0	25.089	43524726	0.0768	0.0766	572235987	100.08
EODET	629141	2	Sigma-Aldrich	08303MEV	0.842	99.4	12.245	32361929	-	-	-	-


Post-Sample IOM												
Container	Standard	CAS	Index #	Manufacturer	Lot #	Density	% purity	Retention time	Area Count	AreaMass	% of Normalized	
Hexane	11543	1	Fisher	13817D 2108257	0.676	99.9	4.340	31911239	0.0652	0.0645	494997450	96.91
Isodecane	14265	3	Sigma-Aldrich	05471R	0.804	99.9	7.615	34912243	0.0652	0.0675	513031456	95.92
Nonane	111842	2	Sigma-Aldrich	ST80788DV	0.724	99.0	14.953	39096996	0.0712	0.0795	154663620	97.38
Decane	124185	2	Sigma-Aldrich	85137H	0.730	99.6	17.321	41354281	0.0729	0.0726	569670320	100.00
Undecane	1120214	3	Sigma-Aldrich	W89029TV	0.740	99.4	19.232	42475211	0.0743	0.0739	514020699	100.93
Dodecane	112403	2	Sigma-Aldrich	W89750AV	0.750	99.0	20.895	44013993	0.0750	0.0742	592799385	104.08
Tridecane	62955	1	Aldrich	W89156RV	0.756	99.0	22.369	44689143	0.0761	0.0761	593209294	104.16
Tetradecane	62954	3	Sigma-Aldrich	ST80371V	0.764	99.0	23.789	44448239	0.0760	0.0762	604011343	100.05
Pentadecane	62953	1	Sigma-Aldrich	101170J	0.769	99.0	25.09	45017745	0.0768	0.0766	602811781	100.08
EODET	629141	2	Sigma-Aldrich	08303MEV	0.842	99.4	12.252	33435150	-	-	-	-

Peripheral Recoveries					FID Retention Times (min)		
Pre-sample	Prepared (g/L)	Measured (g/L)	% Recovery	Raw Counts	49	26	
Pre-sample	EODET	5.814	5.2	104.8	22881929	12.25	12.25
Pre-sample	Triglyme	0.1054	0.1	105.4	548777	31.95	31.95
Post-sample	EODET	5.814	5.4	108.1	33434150		
Post-sample	Triglyme	0.1054	0.1	109.0	5684727		

IOM CHROMATOGRAM

9





10

CONFIGURATION LOOP

- 1) Configure Instrument
- 2) IOM Preparation
- 3) IOM Analysis → Passing QC

11

CONFIGURATION LOOP

- 1) Tweak Configuration
- 2) IOM Preparation
- 3) IOM Analysis → Passing QC

x2

● Failing QC

12

PRE-ANALYSIS CALIBRATION

13

REQUIRED BY EPA

THINGS TO CALIBRATE

A LIST, BY BRAD PARRACK

14

- 1) ???
- 2) ???
- 3) ???

DO YOU KNOW WHAT'S IN YOUR SAMPLE?

15

INDUSTRY LABS



EVERYONE ELSE



EVERYONE ELSE

SAMPLE PRE-SCREENING

- 1) IDENTIFY THE COMPOUNDS IN A SAMPLE
- 2) ESTIMATE THEIR CONCENTRATION TO DETERMINE IF A COMPOUND REQUIRES DIRECT CALIBRATION OR IF A SURROGATE CAN BE USED

ALERTS USER TO PRESENCE OF EXEMPT COMPOUNDS WHICH ARE QUANTIFIED ELSEWHERE AND REMOVED FROM THE VOC COATING CALCULATION


16

SAMPLE PRE-SCREENING

- 1) ~ 3 G OF SAMPLE DILUTED INTO ~ 10 ML OF SOLVENT
- 2) 150 UL OF EGDE (IS) ADDED
- 3) MIXTURE BROUGHT TO 25 ML IN FLASK


17

SCREENING CHROMATOGRAM



18

DEFAULT COMPOUND



TRIGLYME

RESPONSE FACTOR USED TO QUANTIFY PEAKS SEEN:

- 1) DURING SCREENING
- 2) DURING QUANTITATIVE ANALYSIS FOR PEAKS <1 G/L

19

THINGS TO CALIBRATE

A LIST, BY BRAD PARRACK

- 1) TRIGLYME
- 2) HEXANE - PENTADECANE
- 3) ???

20

SAMPLE SPIKING QC

Extraction Markers

IPA

DIURA

Triglyme

Heptane

**DETERMINE
EXTRACTION
EFFICIENCY IN
SOLVENT**

Calibrate

21

THINGS TO CALIBRATE

A LIST, BY BRAD PARRACK

22

- 1) TRIGLYME
- 2) HEXANE - PENTADECANE
- 3) IPA, HEPTANE, DIIBA

OVERHEARD IN THE LAB

23

CALIBRATION TIME!

An Excited Analyst

CALIBRATIONS CURVES

24

4 Point Calibration	2 Point Calibration
Surrogate Spikes	Most Compounds
Difficult Compounds	

CALIBRATIONS CURVES

25

4 Point Calibration	2 Point Calibration
0.1 g/L 1 g/L 5 g/L 15 g/L Forced Through Zero	0.1 g/L 15 g/L Forced Through Zero

4 Point Calibration	2 Point Calibration
---------------------	---------------------

26

Inj. #	Sample	Inj. #	Sample
1	Method Blank	1	Method Blank
2	IOM	2	IOM
3	Method Blank	3	Method Blank
4	15 g/L	4	15 g/L
5	Method Blank	5	Method Blank
6	15 g/L	6	15 g/L
7	Method Blank	7	Method Blank
8	5 g/L	8	0.1 g/L
9	Method Blank	9	Method Blank
10	5 g/L	10	0.1 g/L
11	Method Blank	11	Method Blank
12	1 g/L	12	IOM
13	Method Blank		
14	1 g/L		
15	Method Blank		
16	0.1 g/L		
17	Method Blank		
18	0.1 g/L		
19	Method Blank		
20	IOM		

CALIBRATION CURVES

27

2 POINT CURVE PREP - TRIGLYME & DIIBA

15 g/L Standard			
	Addition	Mass (g)	g/L
25 ml Volumetric Flask (Class A)	-	46.933	-
THF	~10 mL	61.6544	-
EEOE	150 uL	61.7815	5.0536
From stock	-	-	-
DIIBA	397 uL	62.1671	15.2698
Triglyme	397 uL	62.5471	15.0024
QS with THF	To 25 mL	68.921	-

0.1 g/L Standard			
Diluted from 5 g/L standard			
	Addition	Mass (g)	g/L
25 ml Volumetric Flask (Class A)	-	46.9633	-
THF	~10 mL	59.7951	-
EEOE	150 uL	59.8323	5.0575
From 5 g/L stock	500 uL	60.2856	-
DIIBA	-	-	0.1036
Triglyme	-	-	0.0977
QS with THF	To 25 mL	68.9922	-

I FEEL THE NEED.
THE NEED TO SCREEN.

Tom Cruise

28


SAMPLE SCREENING SEQUENCE

Inj. #	Sample
1	Method Blank
2	IOM
3	Method Blank
4	0.1 g/L Exempts Standard
5	0.1 g/L Exempts Standard
6	Method Blank
7	Example Sample
8	Method Blank
9	Example Sample
10	Method Blank
11	0.1 g/L Exempts Standard
12	0.1 g/L Exempts Standard
13	Method Blank
14	IOM

6 SAMPLE MAX

29

SCREENING CHROMATOGRAM



30

WHEN YOU NEED TO HAVE A RESPONSE FACTOR


31

Concentration (g/L)	Action
3+	Calibration Required
1 - 3	Calibration Required -or- Calibrate for Surrogate Required
0.1 - 1	Do Not Calibrate (quant as default compound)
< 0.1	Not counted in sample

Calibrate

SCREENING CHROMATOGRAM

32



SCREENING RESULTS AS TRIGLYME

33

-15 D	FID RT	CMPD Seen	~[g/L]
18	16.83	2-Propanol, 1-butoxy-	5.0
23	18.79	1-Hexanol, 2-ethyl-	0.4
26	21.67	1,3-Pentanediol, 2,2,4-trimet	0.7
27	24.25	Texanol	0.8
28	24.48	Texanol	1.2

REQUIRED CALIBRANT
REQUIRED CALIBRANT
-OR-
REQUIRED SURROGATE
QUANT AS DEFAULT

A NEW LIST OF THINGS TO CALIBRATE 34

A SECOND LIST, BY BRAD PARRACK

- 1) PGBE
- 2) TEXANOL

HAPPY CALIBRATING!

OVERHEARD IN THE LAB


35

ANALYSIS TIME!

An Excited Analyst

SAMPLE ANALYSIS 36

- 1) ~ 30 G OF SAMPLE SPIKED WITH 1 G/L OF EACH EXTRACTION MARKER [IPA/DIIBA/HEPANE/TRIGLYME]
- 2) ~ 3 G OF SAMPLE DILUTED INTO ~ 10 ML OF SOLVENT
- 3) 150 UL OF EGDE (IS) ADDED
- 4) MIXTURE BROUGHT TO 25 ML IN FLASK



37

QUALITY CONTROLS STANDARDS

Continuing Spike Verification CSV	Continuing Cal Verification CCV
Monitors Spike Compounds for Cal Drift	Monitors Target Compounds for Cal Drift
EGDE (IS)	EGDE (IS)
Heptane	PGBE
Triglyme	Texanol
IPA	Target Compounds From Other Samples in Sequence
DIIBA	

38

SAMPLE ANALYSIS SEQUENCE

Inj. #	Sample
1	Method Blank
2	IOM
3	Method Blank
4	CSV
5	CSV
6	Method Blank
7	CCV
8	Method Blank
9	CCV
10	Method Blank
11	Example Sample #1
12	Method Blank
13	Example Sample #1
14	Method Blank
15	Example Sample #2
16	Method Blank
17	Example Sample #2
18	Method Blank
19	CCV
20	Method Blank
21	CCV
22	Method Blank
23	CSV
24	CSV
25	Method Blank
26	IOM

3 SAMPLE MAX

39

Required or Flexible?

40

Flexible	Required
Post Column Split, Two Column System, or Two Discrete Instruments	G43 Column Phase
Solvent (MeOH vs THF)	Column Temp Limit > 300 °C
Consumables	Compound Calibration Prior To Analysis
Split Vs. Splitless	
Column Dimensions	

EASY AS 1,2,3 ... 4,5,6,7,8,9

41

- 1) CONFIGURE INSTRUMENT
- 2) PREP IOM; EVALUATE RESULTS
- 3) RECONFIGURE INSTRUMENT IF NECESSARY
- 4) CALIBRATE DEFAULT COMPOUNDS AND SPIKES
- 5) SCREEN SAMPLES IF UNKNOWN
- 6) CALIBRATE SCREENED/KNOWN COMPOUNDS AND SURROGATES FROM 1+ G/L (CALCULATED AS TRIGLYME)
- 7) PREPARE CCV, CSV; SPIKE/DILUTE SAMPLE
- 8) ANALYZE ON INSTRUMENT
- 9) PROCESS RESULTS

ME, RIGHT NOW

DEMONSTRATION TIME!


Brad Parrack

42

TOUR GROUPS

43

Brad	Tereso	Hanna
David Darling	Stephen Foster	Kristy Rodriguez
Dave Nevisson	Gina Johnson	Kimberly Gutierrez
Chrissy Ford	Brian Morehouse	Jerry Powers
Cathy Willis	Chris Pollack	Chris Nardi
Cidnie Hoang	Barry Cupp	Barry Marcks



44

DATA PROCESSING

DEVELOPING A CALIBRATION LIBRARY

45

CALCULATE RESPONSE FACTORS FOR:

- Internal Standard
- Target Compounds

EACH SOLVENT REQUIRES ITS OWN SET OF CALIBRANTS

4 Point Calibration

2 Point Calibration

46

Inj. #	Sample
1	Method Blank
2	IOM
3	Method Blank
4	15 g/L
5	Method Blank
6	15 g/L
7	Method Blank
8	5 g/L
9	Method Blank
10	5 g/L
11	Method Blank
12	1 g/L
13	Method Blank
14	1 g/L
15	Method Blank
16	0.1 g/L
17	Method Blank
18	0.1 g/L
19	Method Blank
20	IOM

Inj. #	Sample
1	Method Blank
2	IOM
3	Method Blank
4	15 g/L
5	Method Blank
6	15 g/L
7	Method Blank
8	0.1 g/L
9	Method Blank
10	0.1 g/L
11	Method Blank
12	IOM

INTERNAL STANDARD RESPONSE

47

Retention Time	Area Counts	IS Prepared (g/L)
12.272	366039509	5.053
12.271	367831644	5.053
12.27	368661304	5.057
12.272	370512860	5.057

Prepared concentrations and area counts compiled for internal standard

Response Factor for Internal Standard:

$$\frac{\sum \text{Prep'd Concentration} * \text{IS Area}}{\sum \text{Prep'd Concentration}^2}$$

Values are tabulated for Concentrations and Area Counts

X*Y	X^2
1849779195	25.53782
1858835742	25.53782
1864494222	25.57802
1873858415	25.57802
Sum	7446967574 102.2317

INTERNAL STANDARD RESPONSE

48

Response Factor for Internal Standard:

$$\frac{7446967574}{102.2317} = 72843931$$

Response factor applied to area counts to determine percent recoveries of each injection.

Retention Time	Area Counts	Measured g/L	Prepared (g/L)	Percent Recovery
12.272	366039509	5.02	5.053	99.4
12.271	367831644	5.05	5.053	99.9
12.27	368661304	5.06	5.057	100.1
12.272	370512860	5.09	5.057	100.6

TARGET COMPOUND RELATIVE RESPONSE 49

Prepared concentrations and area counts compiled for internal standard and target compound

Prepared (g/L)	Retention Time	Area Counts	IS Prepared (g/L)	IS Retention Time (min)	IS Area Counts
15.270	8.872	1563769174	5.053	12.272	366039509
15.270	8.872	1572094144	5.053	12.271	367831644
0.104	8.755	10723977	5.057	12.27	368661304
0.104	8.757	10758046	5.057	12.272	370512860

Relative Response Factor (RRF) for Target Compound: $\frac{\sum ([\text{Target }]/[\text{IS}]) \times (\text{Target AC}/\text{IS AC})}{\sum ([\text{Target }]/[\text{IS}])}$

TARGET COMPOUND RELATIVE RESPONSE 50

Relative Response Factor for Target Compound:

$$\frac{25.82423}{18.26125} = 1.4142$$

X*Y	X*2
12.90877	9.130206
12.91427	9.130206
0.000596	0.000419
0.000595	0.000419
Sum	25.82423 18.26125

RRF applied to area counts to determine percent recoveries of each injection.

Retention Time	Area Counts	Measured g/L	Prepared (g/L)	Percent Recovery
8.872	1563769174	15.27	15.270	100.0
8.872	1572094144	15.27	15.270	100.0
8.755	10723977	0.10	0.104	100.4
8.757	10758046	0.10	0.104	100.3

RELATIVE RESPONSE FACTORS 51

COMPOUND RRF CAL DATE

1) Triglyme	0.7579	03/15/16
2) DIIBA	1.3850	03/15/16
3) Heptane	2.0510	03/15/16
4) IPA	1.0586	04/06/16
5) Nonane	2.0814	04/06/16
6) Octanol	1.7614	05/17/16
7) 4-Heptanone	1.5783	05/17/16
8) Ethylene Glycol	0.5928	05/17/16
9) Propylene Glycol	0.8091	05/17/16
10) Dipropylene Glycol	0.8981	05/17/16

DETECTION LIMITS

52

Required for all compounds with a RRF lower than Triglyme's RRF

RELATIVE RESPONSE FACTORS

53

COMPOUND	RRF	LOD	DATE
1) Triglyme	0.7579	-	03/15/16
2) DIIBA	1.3850	-	03/15/16
3) Heptane	2.0510	-	03/15/16
4) IPA	1.0586	-	04/06/16
5) Nonane	2.0814	-	04/06/16
6) Octanol	1.7614	-	05/17/16
7) 4-Heptanone	1.5783	-	05/17/16
8) Ethylene Glycol	0.5928	0.004 G/L	05/17/16
9) Propylene Glycol	0.8091	-	05/17/16
10) Dipropylene Glycol	0.8981	-	05/17/16

QUALITY CONTROLS

54

- MB** METHOD BLANK
- IOM** INSTRUMENT OPTIMIZATION MIX
- CSV** CONTINUING SPIKE VERIFICATION
- CCV** CONTINUING CALIBRATION VERIFICATION

CONTINUING SPIKE VERIFICATION 55

CSV
Monitors Spike Compounds for Cal Drift
EGDE (IS)
Heptane
Triglyme
IPA
DIIBA

SPIKES PREPARED AT 1 G/L
PREPARED GRAVIMETRICALLY

CONTINUING SPIKE VERIFICATION 56

Step	Name	grams	net (g)	%Purity	g/L
Tare	Flask	47.0628			
solvent	THF	59.5599	12.4971	99.0	494.8852
CMPD 1	EOEOE	59.6893	0.1294	99.4	5.1449
CMPD 2	Triglyme	59.7148	0.0255	98.7	1.0067
CMPD 3	IPA	59.7386	0.0238	99.99	0.9519
CMPD 4	DIIBA	59.7658	0.0272	99.2	1.0793
CMPD 5	Heptane	59.7870	0.0212	99.0	0.8395
CMPD 6					
CMPD 7					
solvent	THF	68.9512	9.1642	99.0	362.9023

Preparation

CONTINUING SPIKE VERIFICATION 57

peak #	R.T. min	Start min	End min	PK TY	peak height	corr. area	corr. % max.	% of total
1	1.285	1.156	1.330	BV	8743	179832	0.00%	0.000%
2	1.369	1.330	1.596	VB	27721	334469	0.00%	0.000%
3	3.972	3.856	4.159	BB	32991	1196474	0.00%	0.002%
4	4.361	4.283	4.716	BB	18545689	372528076	0.52%	0.470%
5	5.256	5.143	5.436	BV	34846	1786854	0.00%	0.002%
6	6.103	5.436	6.654	PV	622582593	73322864611	100.00%	92.492%
7	6.700	6.654	6.983	VB	6205	585744	0.00%	0.001%
8	7.161	6.996	7.390	BB	21757	1853989	0.00%	0.001%
9	7.539	7.483	7.996	BB	112462008	393568387	0.54%	0.496%
10	8.183	8.090	8.383	BB	15887	577742	0.00%	0.001%
11	8.530	8.410	8.743	BB	38743	2303056	0.00%	0.003%
12	9.628	9.536	9.916	BB	58490	2274521	0.00%	0.003%
13	10.420	10.303	10.650	BB	18425	871900	0.00%	0.001%
14	11.536	11.493	11.616	PV	10553387	429166155	0.99%	0.541%
15	11.751	11.616	12.136	VV	2109407	161329869	0.22%	0.204%
16	12.268	12.136	12.563	VV	9821362	363624514	0.50%	0.459%
17	12.635	12.563	12.870	VB	11579	713796	0.00%	0.001%
18	13.961	13.870	13.907	BV	9911	511855	0.00%	0.000%
19	13.190	13.007	13.596	VV	3826632	224410871	0.31%	0.283%
20	13.655	13.596	13.730	VV	16678	508239	0.00%	0.001%
21	14.963	14.829	15.144	PV	15462512	421896327	0.58%	0.532%
22	15.415	15.228	15.482	VV	9546	626386	0.00%	0.001%
23	15.538	15.482	15.670	VV	7208	476384	0.00%	0.001%
24	15.718	15.670	15.888	VV	7232	278856	0.00%	0.000%
25	15.949	15.888	16.110	VV	13328	520755	0.00%	0.001%

Results From Instrument

CONTINUING SPIKE VERIFICATION 58

RT	Area	Confirmation	Vlookup RF	final VOC g/L
1.19	221335	default	0.758	0.00
1.35	461406	default	0.758	0.00
1.77	151344	default	0.758	0.00
3.36	64809430	default	0.758	0.01
6.06	65579808737	solvent	0.000	0.00
6.86	742690	default	0.758	0.00
7.14	974715	default	0.758	0.00
7.51	12177562	default	0.758	0.01
8.16	538005	default	0.758	0.00
8.51	668143	default	0.758	0.00
9.61	2296732	default	0.758	0.00
12.23	343765372	default	0.758	0.04
17.22	487029	default	0.758	0.00
17.60	366312	default	0.758	0.00
19.78	169408	default	0.758	0.00
21.93	5448231	default	0.758	0.01
22.17	247001	default	0.758	0.00
24.60	467552	default	0.758	0.00
25.74	14096929	default	0.758	0.00
26.45	614906	default	0.758	0.00
28.00	103473809	default	0.758	0.01
28.57	335362	default	0.758	0.00

Compiled Peak Data

CONTINUING SPIKE VERIFICATION 59

RT	Area	Confirmation	Vlookup RF	final VOC g/L
1.19	221335	default	0.758	0.00
1.35	461406	default	0.758	0.01
1.77	151344	default	0.758	0.00
3.36	64809430	IPA	1.059	0.92
6.06	65579808737	solvent	0.000	0.00
6.86	742690	default	0.758	0.01
7.14	974715	default	0.758	0.02
7.51	12177562	Heptane	2.051	0.89
8.16	538005	default	0.758	0.01
8.51	668143	default	0.758	0.01
9.61	2296732	default	0.758	0.05
12.23	343765372	EGDE	6.175	5.57
17.22	487029	default	0.758	0.01
17.60	366312	default	0.758	0.01
19.78	169408	default	0.758	0.00
21.93	5448231	Triglyme	0.758	1.07
22.17	247001	default	0.758	0.00
24.60	467552	default	0.758	0.01
25.74	14096929	default	0.758	0.28
26.45	614906	default	0.758	0.01
28.00	103473809	DIIBA	1.400	1.11
28.57	335362	default	0.758	0.01

Apply RRFs to Appropriate Peaks

CONTINUING SPIKE VERIFICATION 60

CSV Results			
Compound	Prepared	Measured	% Recovery
EGDE (IS)	5.14	5.57	108
Heptane	0.84	0.89	106
Triglyme	1.00	1.07	107
IPA	0.95	0.92	97
DIIBA	1.08	1.11	103

QC Allowance: 90 - 110%

CONTINUING CALIBRATION VERIFICATION 61

CCV should possess target compounds from all samples in the sequence

Same data processing concept as CSV

QC Allowance: 85 - 115%

SAMPLE PROCESSING 62

Step	Name	CAS (required for CMPDs)	weight (g)	net (g)	% Purity	Final wt fx
tare	vial		28.9927			
sample (~35 g)	Sample		57.1552	28.1625	100	0.956
CMPD 1	TRIGLYME	112492	57.4048	0.2496	98.7	0.008
CMPD 2	IPA	67630	57.6621	0.2573	99.99	0.009
CMPD 3	DIIBA	141048	57.9040	0.2419	99.2	0.008
CMPD 4	HEPTANE	142825	58.1513	0.2473	99.0	0.008

Sample Spiking


SAMPLE PROCESSING 63

Step	Name	weight (g)	net (g)	% Purity	Final wt fx	g/L
tare	25 ml flask	51.0738				
solvent	THF	61.6716	10.5978	99.0		
I.S.	EGDE	61.7981	0.1265	99.4	0.006	5.030
spiked sample		65.4544	3.6563	100.0	0.156	141.256
solvent	THF	73.6818	8.2274	99.0		
					Final wt fx	g/L
	TRIGLYME				0.0014	1.236
	IPA				0.0014	1.290
	DIIBA				0.0013	1.204
	HEPTANE				0.0014	1.228

Sample Dilution

SAMPLE PROCESSING

64



SAMPLE PROCESSING

65

peak #	R.T. min	Start min	End min	PK TV	peak height	corr. area	corr. % max	% of total
1	1.205	1.156	1.330	BV	8763	178832	0.00%	0.000%
2	1.369	1.330	1.596	VB	27721	334469	0.00%	0.000%
3	3.972	3.856	4.150	BB	32991	1196474	0.00%	0.002%
4	4.361	4.283	4.716	BB	10945689	37928076	0.52%	0.473%
5	5.256	5.143	5.436	BV	34846	1786854	0.00%	0.002%
6	6.183	5.436	6.654	PV	622582593	73322864611	100.00%	92.492%
7	6.780	6.654	6.983	VB	6205	585744	0.00%	0.001%
8	7.161	6.996	7.390	BB	21757	1853889	0.00%	0.001%
9	7.539	7.483	7.996	BB	11244280	393560387	0.54%	0.496%
10	8.183	8.890	8.383	BB	15887	577742	0.00%	0.001%
11	8.530	8.410	8.743	BB	38743	2383856	0.00%	0.003%
12	9.628	9.536	9.916	BB	58490	2274521	0.00%	0.003%
13	10.420	10.383	10.650	BB	18425	671900	0.00%	0.001%
14	11.536	11.483	11.616	PV	18553367	429166155	0.59%	0.541%
15	11.751	11.616	12.136	VV	2189487	161329869	0.22%	0.204%
16	12.268	12.136	12.563	VV	9821362	363624514	0.58%	0.459%
17	12.635	12.563	12.870	VV	11579	713796	0.00%	0.001%
18	12.961	12.870	13.087	BV	9911	311885	0.00%	0.000%
19	13.190	13.007	13.596	VV	3826632	224418871	0.31%	0.283%
20	13.655	13.596	13.730	VV	16678	580239	0.00%	0.001%
21	14.963	14.829	15.144	PV	15462512	421896327	0.58%	0.532%
22	15.415	15.228	15.482	VV	9546	626386	0.00%	0.001%
23	15.538	15.482	15.670	VV	7288	476384	0.00%	0.001%
24	15.718	15.670	15.888	VV	7232	278856	0.00%	0.000%
25	15.949	15.888	16.110	VV	13328	528755	0.00%	0.001%

Sample Data



SAMPLE PROCESSING

66

Peak No	Spl RT	Spl Area	Confirmation	Vlookup Rf
1	1.77	201680	default	0.758
2	1.92	429948	default	0.758
3	3.24	17363809	default	0.758
4	3.38	8939639	default	0.758
5	6.05	57107053064	default	0.758
6	6.87	929388	default	0.758
7	7.15	635243	default	0.758
8	7.51	65889693	default	0.758
9	8.17	570912	default	0.758
10	9.64	2128814	default	0.758
11	12.24	328099881	default	0.758
12	14.59	247283	default	0.758
13	15.06	225204	default	0.758
14	16.83	112711623	default	0.758
15	17.28	7581679	default	0.758
16	17.61	716721	default	0.758
17	21.94	64183259	default	0.758
18	22.18	359277	default	0.758
19	24.60	356708	default	0.758
20	25.74	11812778	default	0.758
21	26.45	750768	default	0.758
22	27.11	227350	default	0.758
23	28.00	96717162	default	0.758
24	28.58	313003	default	0.758

PEAK ASSESSMENT

All peaks assigned
Triglyme RRF initially



SAMPLE PROCESSING

67

Peak No	Spl RT	Spl Area	Confirmation	Vlookup	Meas RF	Spl VOC g/L	Spl VOC g/L
1	1.77	201680	default	0.758	0.004	0.03	
2	1.92	429948	default	0.758	0.009	0.06	
3	3.24	17363809	IPA	1.059	0.251	1.77	
4	3.38	69339639	IPA	1.059	1.001	7.07	
5	6.05	57107853064	default	0.758	0.000	0.00	
6	6.87	929388	default	0.758	0.019	0.13	
7	7.15	635243	default	0.758	0.013	0.09	
8	7.51	167816779	Heptane	2.051	1.251	8.85	
9	8.17	570912	default	0.758	0.012	0.08	
10	9.64	2128814	default	0.758	0.043	0.30	
11	12.24	329099881	EGDE	6.17E+7	5.334	37.69	
12	14.59	247283	default	0.758	0.005	0.04	
13	15.06	225204	default	0.758	0.005	0.03	
14	16.83	172711923	default	0.758	3.483	24.61	
15	17.26	7591679	default	0.758	0.153	1.08	
16	17.61	716721	default	0.758	0.014	0.10	
17	21.94	64183259	Triglyme	0.758	1.294	9.14	
18	22.18	359277	default	0.758	0.007	0.05	
19	24.60	356708	default	0.758	0.007	0.05	
20	25.74	11812778	default	0.758	0.238	1.68	
21	26.45	750768	default	0.758	0.015	0.11	
22	27.11	227350	default	0.758	0.005	0.03	
23	28.00	96717162	DIIBA	1.400	1.056	7.46	
24	28.58	313033	default	0.758	0.006	0.04	

SPIKE ASSIGNMENT

IS ASSIGNMENT

SAMPLE PROCESSING

68

Sample Spike Results

Compound	Prepared	Measured	% Recovery
EGDE (IS)	5.03	5.33	106
Heptane	1.23	1.25	102
Triglyme	1.24	1.29	104
IPA	1.29	1.25	97
DIIBA	1.20	1.06	88

QC Allowance: 85 - 115%

SAMPLE PROCESSING

69

Pre Blank RT	Pre Blank AC	Pre CSV RT	Pre CSV AC	Peak No	Spl RT	Spl Area	Confirmation	Vlookup	Meas RF	Spl VOC g/L	Spl VOC g/L
1.77	176241	1.77	194712	1	1.77	201680	default	0.758	0.004	0.03	
		3.22	17868001	3	3.24	17363809	IPA	1.059	0.251	1.77	
		3.37	48640561	4	3.38	69339639	IPA	1.059	1.001	7.07	
6.06	67951221220	6.07	67977122016	5	6.05	57107853064	default	0.758	0.000	0.00	
		6.86	764954	6	6.87	929388	default	0.758	0.019	0.13	
7.15	712058	7.14	1096999	7	7.15	635243	default	0.758	0.013	0.09	
		7.51	129162718	8	7.51	55939593	Heptane	2.051	0.417	2.95	
8.17	552248	8.16	560389	9	8.17	570912	default	0.758	0.012	0.08	
9.61	2331069	9.61	2327567	10	9.64	2128814	default	0.758	0.043	0.30	
12.24	353833505	12.23	353199810	11	12.24	329099881	EGDE	6.17E+7	5.334	37.69	
				12	14.59	247283	default	0.758	0.005	0.04	
				13	15.06	225204	default	0.758	0.005	0.03	
				14	16.83	172711923	default	0.758	3.483	24.61	
				15	17.26	7591679	default	0.758	0.153	1.08	
17.60	1145014	17.60	1170358	16	17.61	716721	default	0.758	0.014	0.10	
		21.94	5606829	17	21.94	64183259	Triglyme	0.758	1.294	9.14	
		22.17	245218	18	22.18	359277	default	0.758	0.007	0.05	
24.60		24.60	478330	19	24.60	356708	default	0.758	0.007	0.05	
25.74	14677563	25.74	14755526	20	25.74	11812778	default	0.758	0.238	1.68	
		26.45	721685	21	26.45	750768	default	0.758	0.015	0.11	
		28.00	105953275	22	27.11	227350	default	0.758	0.005	0.03	
		28.58	338987	23	28.00	96717162	DIIBA	1.400	1.056	7.46	
				24	28.58	313033	default	0.758	0.006	0.04	

SAMPLE PROCESSING

70

Peak No	Spl RT	Spl Area	Confirmation	Vlookup RF	Spl VOC g/L
1	1.77	201680	solvent	0.000	0.00
2	1.92	429948	default	0.758	0.06
3	3.24	17363809	IPA	1.059	1.77
4	3.38	89339639	IPA	1.059	7.07
5	6.05	57107853064	solvent	0.000	0.00
6	6.87	929388	default	0.758	0.13
7	7.15	635243	solvent	0.000	0.00
8	7.51	55939593	Heptane	2.051	2.96
9	8.17	570912	solvent	0.000	0.00
10	9.64	2128814	solvent	0.000	0.00
11	12.24	329099881	EGDE	6.17E+7	37.69
12	14.59	247283	default	0.758	0.04
13	15.06	225204	default	0.758	0.03
14	16.83	172711923	default	0.758	24.61
15	17.26	7591679	default	0.758	1.08
16	17.61	716721	solvent	0.000	0.00
17	21.94	84183259	Triglyme	0.758	9.14
18	22.18	359277	default	0.758	0.05
19	24.60	356708	default	0.758	0.05
20	25.74	11812778	solvent	0.000	0.00
21	26.45	750768	default	0.758	0.11
22	27.11	227350	default	0.758	0.03
23	28.00	96717162	DIBA	1.400	7.46
24	28.58	313033	default	0.758	0.04

SAMPLE PROCESSING

71

Results				
Compound Seen	Concentration as Triglyme	RRF Used	RRF	Final Concentration
PGBE	24.61	PGBE	1.179	15.82
2-ethyl hexanol	1.08	Octanol	0.84	0.48

SAMPLE PROCESSING

72

Peak No	Spl RT	Spl Area	Confirmation	Vlookup RF	Spl VOC g/L
1	1.77	201680	solvent	0.000	0.00
2	1.92	429948	default	0.758	0.06
3	3.24	17363809	IPA	1.059	1.77
4	3.38	89339639	IPA	1.059	7.07
5	6.05	57107853064	solvent	0.000	0.00
6	6.87	929388	default	0.758	0.13
7	7.15	635243	solvent	0.000	0.00
8	7.51	55939593	Heptane	2.051	2.96
9	8.17	570912	solvent	0.000	0.00
10	9.64	2128814	solvent	0.000	0.00
11	12.24	329099881	EGDE	6.17E+7	37.69
12	14.59	247283	default	0.758	0.04
13	15.06	225204	default	0.758	0.03
14	16.83	172711923	PGBE	1.179	15.82
15	17.26	7591679	Octanol	1.716	0.48
16	17.61	716721	solvent	0.000	0.00
17	21.94	84183259	Triglyme	0.758	9.14
18	22.18	359277	default	0.758	0.05
19	24.60	356708	default	0.758	0.05
20	25.74	11812778	solvent	0.000	0.00
21	26.45	750768	default	0.758	0.11
22	27.11	227350	default	0.758	0.03
23	28.00	96717162	DIBA	1.400	7.46
24	28.58	313033	default	0.758	0.04

SAMPLE PROCESSING

73

Final Sample Results (g/L)	
Toal VOC	82.62
Surrogate Spikes	66.08
Solvent (not counted)	0.32
Total Left As Default	0.24
Net VOC Material	16.54

SAMPLE PROCESSING

74

Final Sample Results	
Total Left As Default	0.24

Limit: 5 g/L or 10%; whichever is larger

SAMPLE PROCESSING SUMMARY

75

- Spike Samples with Extraction Markers
- Dilute Sample
- Inject Sample
- Import Data, Check Integrations, Assign Every Peak the Triglyme RRF
- Assign RF to EOEOE peak; RRFs to Extraction Marker peaks
- Remove Solvent Peaks from Summation
- Identify Peaks > 1 g/L; Assign Direct RRF for Peaks > 3 g/L
- Assign Direct RRF for Peaks > 1 and Less than 3 g/L -or- Assign Acceptable Surrogate RRF
- Report Total VOC Material and QC for IOMs, CCVs, and CSVs

SURROGATE RULES

76


Between 1 - 3 g/L surrogate compounds may be utilized to quantify so long as they are within the following rules:

- 1) Surrogate compounds possess the same functional groups
- 2) Surrogate compounds are within 1 Carbon


SURROGATE RULES

77

Example:
2-ethyl hexanol [2 g/L]




Acceptable Surrogate
Octanol




SURROGATE RULES

78

Example:
Toluene [1.5 g/L]




Acceptable Surrogate:
Ethyl Benzene




SURROGATE RULES

79

Example:
Butyl Acetate [0.75 g/L]



Acceptable Surrogate:
Triglyme
(< 1 g/L)



ERROR BAND TREATMENT

80

Designed to limit error to
5 g/L VOC Material

ERROR BAND

81

EXAMPLE COATING

VOC Material:	20.8 g/L
Wt % Water:	45.6
NV %:	52.7
Density:	1.2 g/mL
VOC Coating:	48.5 g/L

82

VOC Material: 20.8 g/L


VOC Coating: 48.5 g/L

VOC Material: 15.8 g/L
(± 5 g/L) 25.8 g/L

**VOC Coating Range
37 - 60 g/L**



83



84
