

## Example Calibration Sets

Set 1 (0.1/1/10/15)	
Compound	RT

EGDE	15.9
IPA	6-7.6
Heptane	11.2
DIIBA	31.8
TRIG	25.3

Set 2 (0.1/15)
----------------

Toluene	15.0
Propylene Glycol	16.8
Nonane	18.5
Decane	20.8
Dipropylene Glycol	22.8/23.3/23.8
Texanol	27.6/27.8
Benzophenone	31.9

Set 3 (0.1/15)
----------------

Triethylamine	10.9
p-Xylene	18.4
2-Heptanone	19.3
Dipropylene Glycol Methyl Ether	21.6/21.7/21.8/21.9
Dodecane	24.3
TXIB	30.4

## Example Calibration Sets

Set 4 (0.1/15)
----------------

2-Butanone Oxime	16.5
Mesitylene	20.7
Undecane	22.7
NMP	23.2
Dipropylene Glycol Butyl Ether	25.5/25.6

Set 5 (0.2/0.5)
-----------------

Tetraethylene Glycol	30.1
Pentaethylene Glycol	35.6